



BIMM 143
Genome Informatics II


Lecture 14

Barry Grant
UC San Diego


<http://thegrantlab.org/bimm143>

RNA sequencing overview

In vivo



In vitro



In silico



Goal: RNA quantification, transcript discovery, variant identification


Mapping/Alignment



Quantification

Absolute read counts	15	5	15	(35)
Normalized read counts	$RPKM = \frac{totalTranscriptReads}{mappedReads(millions) \times transcriptLength(Kb)}$			(0.7)

Transcript discovery




Variant discovery




RNA Sequencing

The absolute basics

Normal Cells




Mutated Cells




- The **mutated cells** behave differently than the **normal cells**
- We want to know what genetic mechanism is causing the difference
- One way to address this is to examine differences in gene expression via RNA sequencing...


Normal Cells




Mutated Cells






Each cell has a bunch of chromosomes



Normal Cells




Mutated Cells





Each chromosome has a bunch of genes


Normal Cells




Mutated Cells






Some genes are active more than others




Normal Cells




Mutated Cells





Normal Cells




Mutated Cells





HTS tells us which genes are active, and how much they are transcribed!




Normal Cells



Mutated Cells




We use RNA-Seq to measure gene expression in normal cells ...




... then use it to measure gene expression in mutated cells



Normal Cells




Mutated Cells




Then we can compare the two cell types to figure out what is different in the mutated cells!



Normal Cells




Mutated Cells



Gene2

Gene3



Differences apparent for Gene 2 and
to a lesser extent Gene 3

3 Main Steps for RNA-Seq:

1) Prepare a sequencing library

(RNA to cDNA conversion via reverse transcription)

2) Sequence

(Using the same technologies as DNA sequencing)

3) Data analysis


(Often the major bottleneck to overall success!)

We will discuss each of these steps - but we will focus on step 3 today!

Our last class got us to the start of step 3!

Gene	WT-1	WT-2	WT-3	...
A1BG	30	5	13	...
AS1	24	10	18	...
...

We **sequenced, aligned, counted** the reads per gene in each sample to arrive at our **count table/matrix**



Inputs

Control

Reads R1
FastQ

Treatment

Reads R1
FastQ

Inputs

Control

Reads R1
FastQ



Reads R2
[optional]
FastQ


Treatment


Reads R1
FastQ


Reads R2
[optional]
FastQ

Optional Replicates










...Now what?

This is where we stopped last day



Do it Yourself!

Install DESeq2

[Bioconductor Setup Link](#)

```
install.packages("BiocManager")  
BiocManager::install()  
  
# For this class, you'll also need DESeq2:  
BiocManager::install("DESeq2")
```

Note: Answer NO to prompts to install from source or update...

Do it Yourself!

```
Old packages: 'devtools', 'dplyr', 'DT', 'ggplot2', 'ggpubr',  
'lattice', 'MASS', 'Matrix', 'mclust', 'mgcv', 'openssl',  
'packrat', 'pkgload', 'ps', 'psych', 'raster', 'rcmdcheck',  
'Rcpp', 'remotes', 'rsconnect', 'sessioninfo', 'shiny',  
'shinythemes', 'survival', 'tidyr', 'tinytex', 'xfun'
```

```
Update all/some/none? [a/s/n]:
```

n

Install DESeq2


[Bioconductor Setup Link](#)

```
install.packages("BiocManager")  
BiocManager::install()  
  
# For this class, you'll also need DESeq2:  
BiocManager::install("DESeq2")
```

Note: Answer **NO** to prompts to install from source or update...

Background to Today's Data

Glucocorticoids inhibit inflammatory processes and are often used to treat **asthma** because of their anti-inflammatory effects on airway smooth muscle (ASM) cells.



Mechanism?

- Data from: Himes *et al.* "[RNA-Seq Transcriptome Profiling Identifies CRISPLD2 as a Glucocorticoid Responsive Gene that Modulates Cytokine Function in Airway Smooth Muscle Cells.](#)" PLoS ONE. 2014 Jun 13;9(6):e99625.

Background to Today's Data

- The anti-inflammatory effects of glucocorticoids on airway smooth muscle (ASM) cells has been known for some time but the underlying molecular mechanisms are unclear.
- Himes *et al.* used **RNA-seq** to profile gene expression changes in 4 ASM cell lines treated with **dexamethasone** (a common synthetic glucocorticoid).
- Used Tophat and Cufflinks and found many differentially expressed genes. Focus on CRISPLD2 that encodes a secreted protein involved in lung development
- SNPs in CRISPLD2 in previous GWAS associated with inhaled corticosteroid resistance and bronchodilator response in asthma patients.
- Confirmed the upregulated CRISPLD2 with qPCR and increased protein expression with Western blotting.

Data pre-processing

- Analyzing RNA-seq data starts with sequencing reads.
- Many different approaches, see references on class website.
- Our workflow (previously done):
 - Reads downloaded from GEO ([GSE:GSE52778](#))
 - Quantify transcript abundance ([kallisto](#)).
 - Summarize to gene-level abundance ([txImport](#))
- Our starting point is a **count matrix**: each cell indicates the number of reads originating from a particular **gene** (in rows) for each **sample** (in columns).

counts + metadata

1

countData

gene	ctrl_1	ctrl_2	exp_1	exp_2
geneA	10	11	56	45
geneB	0	0	128	54
geneC	42	41	59	41
geneD	103	122	1	23
geneE	10	23	14	56
geneF	0	1	2	0
...

2

colData


id	treatment	sex	...
ctrl_1	control	male	...
ctrl_2	control	female	...
exp_1	treated	male	...
exp_2	treated	female	...

colData describes metadata about the *columns* of countData

countData is the count matrix
(Number of reads coming from each
gene for each sample)

N.B. First column of colData must match column names (i.e. **sample names**) of countData (-1st)

Counting is (relatively) easy:



Do it Yourself!

Hands-on time!


https://bioboot.github.io/bimm143_F19/lectures/#14

Fold change (log ratios)

- To a **statistician** fold change is sometimes considered **meaningless**. Fold change can be large (e.g. >>two-fold up- or down-regulation) without being statistically significant (e.g. based on probability values from a t-test or ANOVA).
- To a **biologist** fold change is almost always considered **important** for two reasons. First, a very small but statistically significant fold change might not be relevant to a cell's function. Second, it is of interest to know which genes are most dramatically regulated, as these are often thought to reflect changes in biologically meaningful transcripts and/or pathways.

Volcano plot

A common summary figure used to highlight genes that are both significantly regulated and display a high fold change



A volcano plot shows fold change (x-axis) versus $-\log$ of the p -value (y-axis) for a given transcript. The more significant the p -value, the larger the $-\log$ of that value will be. Therefore we often focus on 'higher up' points.


```
# Setup your point color vector
```

```
> mycols <- rep("gray", nrow(res01))
```

```
> mycols[ abs(res01$log2FoldChange) > 2 ] <- "red"
```

```
> inds <- (res01$padj < 0.01) & (abs(res01$log2FoldChange) > 2 )
```


```
> mycols[ inds ] <- "blue"
```

```
# Volcano plot with custom colors
```

```
> plot( res01$log2FoldChange, -log(res01$padj),  
       col=mycols, ylab="-Log(P-value)",  
       xlab="Log2(FoldChange)" )
```


```
> abline(v=c(-2,2), col="gray", lty=2)
```

```
> abline(h=-log(0.1), col="gray", lty=2)
```



Volcano Plot

Fold change vs P-value



Recent developments in RNA-Seq

- **Long read sequences:**
 - PacBio and Oxford Nanopore [[Recent Paper](#)]
- **Single-cell RNA-Seq:** [[Review article](#)]
 - Observe heterogeneity of cell populations
 - Detect sub-population
- **Alignment-free quantification:**
 - Kallisto [[Software link](#)]
 - Salmon [[Software link](#), [Blog post](#)]

Additional Reference Slides


Public RNA-Seq data sources

- **Gene Expression Omnibus (GEO):**
 - <http://www.ncbi.nlm.nih.gov/geo/>
 - Both microarray and sequencing data
- **Sequence Read Archive (SRA):**
 - <http://www.ncbi.nlm.nih.gov/sra>
 - All sequencing data (not necessarily RNA-Seq)
- **ArrayExpress:**
 - <https://www.ebi.ac.uk/arrayexpress/>
 - European version of GEO
- All of these have links between them

Paired-End Reads



Alignment to the Reference Sequence



Paired-end sequencing enables both ends of the DNA fragment to be sequenced. Because the distance between each paired read is known, alignment algorithms can use this information to map the reads over repetitive regions more precisely. This results in much better alignment of the reads, especially across difficult-to-sequence, repetitive regions of the genome.

Taken From: <https://www.illumina.com/science/technology/next-generation-sequencing/paired-end-vs-single-read-sequencing.html>

Count Normalization

- Normalization is required to make comparisons in gene expression
 - Between 2+ genes in one sample
 - Between genes in 2+ samples
- Genes will have more reads mapped in a sample with high coverage than one with low coverage
 - $2x$ depth \approx $2x$ expression
- Longer genes will have more reads mapped than shorter genes
 - $2x$ length \approx $2x$ more reads

Normalization: RPKM, FPKM & TPM

- **N.B.** Some tools for differential expression analysis such as edgeR and DESeq2 want raw read counts - i.e. non normalized input!
- However, often for your manuscripts and reports you will want to report normalized counts
- RPKM, FPKM and TPM all aim to normalize for sequencing depth and gene length. For the former:
 - Count up the total reads in a sample and divide that number by 1,000,000 - this is our “per million” scaling.
 - Divide the read counts by the “per million” scaling factor. This normalizes for sequencing depth, giving you reads per million (RPM)
 - Divide the RPM values by the length of the gene, in kilobases. This gives you RPKM.

- FPKM was made for paired-end RNA-seq
- With paired-end RNA-seq, two reads can correspond to a single fragment
- The only difference between RPKM and FPKM is that FPKM takes into account that two reads can map to one fragment (and so it doesn't count this fragment twice).

- **TPM** is very similar to RPKM and FPKM. The only difference is the order of operations:
 - First divide the read counts by the length of each gene in kilobases. This gives you reads per kilobase (RPK).
 - Count up all the RPK values in a sample and divide this number by 1,000,000. This is your “per million” scaling factor.
 - Divide the RPK values by the “per million” scaling factor. This gives you TPM.
- Note, the only difference is that you normalize for gene length first, and then normalize for sequencing depth second.

- When you use TPM, the sum of all TPMs in each sample are the same.
- This makes it easier to compare the proportion of reads that mapped to a gene in each sample.
- In contrast, with RPKM and FPKM, the sum of the normalized reads in each sample may be different, and this makes it harder to compare samples directly.